

# resDetect<sup>TM</sup> Human Interferon-γ (IFN-γ) ELISA Kit (Residue Testing)

(Enzyme-Linked Immunosorbent Assay)

Catalog Number: RES-A028

Pack Size: 96 tests

IMPORTANT: Please carefully read this manual before performing your experiment.

For Research Use Only. Not For Use in Diagnostic or Therapeutic Procedures



## **INTENDED USE**

The Human Interferon- $\gamma$  (IFN- $\gamma$ ) ELISA Kit (Residue Testing) was developed for the detection and quantitative determination of GMP human IFN- $\gamma$  in samples from CAR-T product preparation processing. It is intended for research use only (RUO).

## **BACKGROUND**

Interferon-gamma (IFN-γ/IFNG) is a dimerized soluble cytokine that is the only member of the type II class of interferon. This interferon was originally called macrophage-activating factor, a term now used to describe a larger family of proteins to which IFN-γ belongs. IFN-γ has been used in a wide variety of clinical indications. Interferon-gamma (IFN-γ) is a central regulator of the immune response and signals via the Janus Activated Kinase (JAK)-Signal Transducer and Activator of Transcription (STAT) pathway. Interferon gamma has broader roles in activation of innate and adaptive immune responses to viruses and tumors, in part through upregulating transcription of genes involved in cell cycle regulation, apoptosis, and antigen processing/presentation. Despite this, rodent and human trophoblast cells show dampened responses to IFNG that reflect the resistance of these cells to IFNG-mediated activation of major histocompatibility complex (MHC) class II transplantation antigen expression.

To support the development of CAR-T drugs, ACROBiosystems independently developed human Interferon- $\gamma$  (IFN- $\gamma$ ) ELISA Residue Testing kit via rigorous methodological validation, which is used for detection of GMP human IFN- $\gamma$  in samples from CAR-T product preparation processing for evaluation the quality of CAR-T products in drug development and CMC quality control stages.

# PRINCIPLE OF THE ASSAY

This assay kit is used to measure the levels of human Interferon- $\gamma$  (IFN- $\gamma$ ) by employing a standard sandwich-ELISA format. The micro-plate in the kit has been pre-coated with Anti-IFN- $\gamma$  Antibody. Firstly, add the standard samples provided in kit and your samples to the plate, incubate and wash the wells. Then add the Biotin-Anti-IFN- $\gamma$  Antibody to the plate and form Antibody-antigen-biotinylated antibody complex, incubate and wash the wells. Next add Streptavidin-HRP to the plate, incubate and wash the wells. At last, load the substrate into the wells and monitor solution color from blue to yellow. The reaction is stopped by the addition of a stop solution and the intensity of the absorbance can be measured at 450nm and 630nm. The OD Value reflects the amount of IFN- $\gamma$  bound.

1 / 12



## **PRECAUTIONS**

- 1. This kit is for research use only and is not for use in diagnostic or therapeutic applications.
- 2. The kit is suitable for cell supernatant, serum and plasma samples.
- 3. Do not use reagents past their expiration date.
- 4. Do not mix or substitute reagents with those from other kits or other lot number kits.
- 5. If samples generate values higher than the highest standard, dilute the samples with the appropriate calibrator diluent and repeat the assay. If cell supernatant samples need step dilution, except for the final dilution with diluent, other intermediate dilutions can be in cell culture medium.
- 6. Differences in test results can be caused by a variety of factors, including laboratory operator, pipette usage, plate washing technique, reaction time or temperature, and kit storage.
- 7. This kit is designed to remove or reduce some endogenous interference factors in biological samples, and not all possible influencing factors have been removed.

## MATERIALS PROVIDED

Table 1. Materials provided

Catalan	Comments	Size	E4	Storage	
Catalog	Components	(96 tests)	Format	Unopened	Opened
RES028-C01	Pre-coated Anti-IFN-γ Antibody Microplate	1 plate	Solid	2-8°C	2-8°C
RES028-C02	Human IFN-γ Standard	20 μg	Power	2-8°C	-70°C
RES028-C03	Biotin-Anti-IFN-γ Antibody	50 μL	Liquid	2-8°C	2-8°C
RES028-C04	Streptavidin-HRP	50 μL	Liquid	2-8°C, avoid light	2-8°C, avoid light
RES028-C05	10xWashing Buffer	50 mL	Liquid	2-8°C	2-8°C
RES028-C06	2xDilution Buffer	50 mL	Liquid	2-8°C	2-8°C
RES028-C07	Substrate Solution	12 mL	Liquid	2-8°C, avoid light	2-8°C, avoid light
RES028-C08	Stop Solution	7 mL	Liquid	2-8°C	2-8°C

2 / 12

US and Canada: Tel: +1 800-810-0816 Web: <a href="http://www.acrobiosystems.com">http://www.acrobiosystems.com</a>
Asia and Pacific: Tel: +86 400-682-2521 E-mail: <a href="mailto:order@acrobiosystems.com">order@acrobiosystems.com</a>



## **SRORAGE**

- 1. The unopened kit should be stored at 2°C -8°C upon receiving. Find the expiration date on the outside packaging and do not use reagents past their expiration date.
- 2. The opened kit should be stored per Table 1. The shelf life is 30 days from the date of opening.

## REAGENTS/EQUIPMENT NEEDED BUT NOT SUPPLIED

Single or multi-channel micropipettes and pipette tips: need to meet  $10 \mu L$ ,  $300 \mu L$ ,  $1000 \mu L$  injection requirements;

37°C Incubator;

Single or dual wavelength microplate reader with 450nm and 630nm filter;

Tubes: 1.5mL,10mL;

Timer;

Reagent bottle;

Deionized or distilled water.

## **REAGENT PREPARATION**

Bring all reagents and samples to room temperature (20°C-25°C) before use. If crystals have formed in buffer solution, place the sample in an 37°C incubator until the crystals have completely dissolved and bring the solution back to room temperature before use.

According to Table 2, prepare the provided lyophilized product into a storage solution with ultrapure water, dissolve at room temperature for 15 to 30 minutes, and mix by gently pipetting, avoiding vigorous shaking or vertexing. The reconstituted storage solution should be stored at -70°C. It is recommended that the number of freezing and thawing should not exceed 1 time, and the size of the aliquot should not be less than  $10 \mu g$ .

**Note:** Considering inevitable minor quantitation variations between protein batches, it is also reasonable to generate the standard curve with specific lot of proteins used for current production for even better accuracy.

**Table 2. Preparation method** 

ID	Components	Size (96 T)	Storage solution concentration.	Reconstituted water Vol.
RES028-C02	Human IFN-γ Standard	20 μg	200 μg/mL	100 μL

3 / 12

US and Canada:

Asia and Pacific:

Tel: +1 800-810-0816

Tel: +86 400-682-2521

E-mail: order@acrobiosystems.com



## RECOMMENDED SAMPLE PREPARATION

#### 1. Working Solution Preparation

1.1 Preparation of 1×Washing Buffer:

Dilute 50 mL 10×Washing Buffer with ultrapure water/deionized water to 500 mL.

1.2 Preparation of 1×Dilution Buffer:

Dilute 50 mL 2×Dilution Buffer with 1×Washing Buffer to 100 mL.

1.3 Preparation of Biotin-Anti-IFN-γ Antibody working fluid:

Dilute Biotin-Anti-IFN-γ Antibody at 1:625 with 1×Dilution Buffer. Please prepare it for one-time use only.

1.4 Preparation of Streptavidin-HRP working fluid:

Dilute Streptavidin-HRP at 1:2000 with 1×Dilution Buffer. The prepared working fluid should avoid light. Please prepare it for one-time use only.


- 1.5 Sample preparation
- a. If the sample to be tested is the serum or plasma, dilute test sample at 1:2 with 1×Dilution Buffer. The volume ratio of sample to diluent is 1:1.
- b. If the sample to be tested is the cell supernatant, dilute test sample at 1:2 with 1×Dilution Buffer. The volume ratio of sample to diluent is 1:1.

### 2. Preparation of Standard curve

The concentration of the reconstituted human IFN-γ Calibrator (RES028-C02) is 200 μg/mL, prepare (Std.-0) by diluting 10 μL the reconstituted human IFN-γ Calibrator into 990 μL Sample Dilution Buffer, mix gently well. Then prepare Std.-1' by diluting 10 μL Std.-0 into 990 μL Sample Dilution Buffer. At last, prepare the highest concentration of standard curve, **Std.-1** (200 pg/mL), by diluting 10 μL Std.-1' into 990 μL Sample Dilution Buffer. Prepare 1:1 serial dilutions for the standard curve as follows: Pipette 300 μL of Sample Dilution Buffer into each tube. Make sure to mix well every time. Sample Dilution Buffer serves as blank.

4 / 12





## 3. Add Samples

Add 100 µL Calibrator and samples to each well. For blank Control wells, please add 100 µL Dilution Buffer.

Note: It is recommended to set doable holes for samples and standard curves to be tested.

#### 4. Incubation

Seal the plate with microplate sealing film and incubate at room temperature for 1 hour.

#### 5. Washing

Remove the remaining solution by aspiration, add 300 µL of 1×Washing Buffer to each well, soak for 10s, remove any remaining 1×Washing Buffer: by aspirating or decanting, invert the plate and blot it against paper towels. Repeat the wash step above for three times.

### 6. Add Biotin-Anti-IFN-γ Antibody

For all wells, add 100  $\mu$ L Biotin-Anti-IFN- $\gamma$  Antibody (dilute at 1:625) working solution. Please prepare it for one-time use only.

#### 7. Incubation

Seal the plate with microplate sealing film and incubate at room temperature for 1 hour.

#### 8. Washing

Repeat step 5.

5 / 12

US and Canada:

Asia and Pacific:

Tel: +1 800-810-0816



## 9. Add Streptavidin-HRP

For all wells, add  $100 \mu L$  Streptavidin-HRP (dilute at 1:2000) working solution. Please prepare it for one-time use only, avoid light.

### 10. Incubation

Seal the plate with microplate sealing film and incubate at room temperature for 30 min.

#### 11. Washing

Repeat step 5.

#### 12. Substrate Reaction

Add 100 µL Substrate Solution to each well. Seal the plate with microplate sealing film and incubate at room temperature for 20 min, avoid light.

#### 13. Termination

Add 50 µL Stop Solution to each well, and tap the plate gently to allow thorough mixing.

*Note:* The color in the wells should change from blue to yellow.

## 14. Data Recording

Read the absorbance at 450 nm and 630 nm using UV/Vis microplate spectrophotometer within 10 minutes.

**Note**: To reduce the background noise, subtract the value read at  $OD_{450nm}$  with the value read at  $OD_{630 nm}$ .

## **CALCULATION OF RESULTS**


- 1. Calculate the mean absorbance for each standard, control and sample and subtract average zero standard optical density (O.D.).
- 2. The standard curve is plotted with the standard concentration as x-axis and the calibrated absorbance value as y-axis. Four parameters logistic are used to draw the standard curve and calculate the sample concentration.
- 3. Normal range of Standard curve:  $R^2 \ge 0.9900$ .
- 4. Detection range: 3.125 pg/mL-200 pg/mL. If the OD value of the sample to be tested is higher than 200 pg/mL, the sample shall be diluted with dilution buffer and assay repeated. If the OD value of the sample to be tested is lower than 3.125pg/mL, the sample should be reported.

6 / 12

US and Canada: Tel: +1 800-810-0816 Web: <a href="http://www.acrobiosystems.com">http://www.acrobiosystems.com</a>
Asia and Pacific: Tel: +86 400-682-2521 E-mail: <a href="mailto:order@acrobiosystems.com">order@acrobiosystems.com</a>



# **OUICK GUILD**






## **TYPICAL DATA**

For each experiment, a standard curve needs to be set for each micro-plate, and the specific OD value may vary depending on different laboratories, testers, or equipments. The following example data is for reference only. The sample concentration was calculated based on the results of the standard curve.

Conc.(pg/mL)	O.D1	O.D2	Average	Corrected
200	2.517	2.488	2.557	2.488
100	1.630	1.628	1.690	1.621
50	0.974	0.989	0.997	0.928
25	0.531	0.534	0.540	0.471
12.5	0.316	0.301	0.314	0.245
6.25	0.181	0.177	0.186	0.117
3.13	0.118	0.114	0.121	0.052
0	0.062	0.063	0.069	/



## **SENSITIVITY**

The minimum detectable concentration of human IFN-γ is 0.382 pg/mL. The minimum detectable concentration was determined by adding twice standard deviations to the mean optical density value of twenty zero standard replicates and calculating the corresponding concentration.

### **PRECISION**

### 1. Intra-assay Precision

Three samples of known concentration were tested twenty times on one plate to assess intra-assay precision.

### 2. Inter-assay Precision

Three samples of known concentration were tested in three separate assays to assess inter-assay precision.

	Intra-assay Precision			Inter-assay Precision		
Sample	1	2	3	1	2	3
n	20	20	20	3	3	3
Mean (pg/mL)	181.751	55.235	9.261	181.539	55.270	9.339
SD	6.884	2.712	0.488	3.038	0.086	0.234
CV (%)	3.8	4.9	5.3	1.7	0.2	2.5

Note: The example data is for reference only.

8 / 12

US and Canada: Tel: +1 800-810-0816 Web: <a href="http://www.acrobiosystems.com">http://www.acrobiosystems.com</a>
Asia and Pacific: Tel: +86 400-682-2521 E-mail: <a href="mailto:order@acrobiosystems.com">order@acrobiosystems.com</a>



# **RECOVERY**

Three parts of blank serum were added with different concentrations of human IFN- $\gamma$ , and the serum without human IFN- $\gamma$  was used as background to calculate the recovery rate. The range of the recovery rate is 87.8-109.7%, and the average recovery is 98.5%.

Sample Type	Average % Recovery	Range
Serum(n=5)	98.5%	87.8-109.7%

## **LINEARITY**

To assess the linearity of the assay, samples spiked with high concentrations of human IFN-γ were serially diluted with calibrator diluent to produce samples with values within the dynamic range of the assay.

		Cell culture medium (DMEM)	Cell culture medium (1640)	Serum
Average Recovery (%)		93.5	97.1	99.9
1:2	Range (%)	89.2-94.3	95.2-99.2	97.5-103.9
1.4	Average Recovery (%)	97.4	98.6	104.2
1:4	Range (%)	94.8-99.9	97.2-100.6	102.3-106.6
1.0	Average Recovery (%)	101.6	99.7	106.4
1:8	Range (%)	99.4-105.1	97.2-101.2	104.5-107.8
1.16	Average Recovery (%)	105.2	102.1	108.9
1:16	Range (%)	102.2-107.6	98.3-104.3	104.5-113.8

*Note*: The example data is for reference only.

### **SPECIFICITY**

This assay recognizes natural and recombinant human IFN-γ. No cross-reactivity was observed when this kit was used to analyze the following recombinant cytokines.

Human					
IL-2	IL-11	GM-CSF	L1R		
IL-3	IL-12B	G-CSF	VEGF165		
IL-4	IL-15	M-CSF	Anti-CD3		
IL-5	IL-17A	SCF	Anti-CD28		
IL-6	IL-18	BMP-2	Anti-CD137		
IL-7	IFN-alpha 1	FGF basic	Thrombopoietin-TPO		
IL-8	TNF-α	Flt-3 Ligand	TGF-beta 1		
IL-10					

9 / 12

 US and Canada:
 Tel: +1 800-810-0816
 Web: <a href="http://www.acrobiosystems.com">http://www.acrobiosystems.com</a>

 Asia and Pacific:
 Tel: +86 400-682-2521
 E-mail: order@acrobiosystems.com



# **NTERFERING SUBSTANCES**

Verify potential matrix effects by adding different levels of DMSO and HSA to the diluted buffer.


Additive	Tolerated concentration
DMSO	10%
HSA	5%

## **CALIBRATION**

This immunoassay is calibrated against a highly purified *E. coli*-expressed recombinant human IFN-γ (87/586). Reference Reagent is calibrated by NIBSC/WHO in April 2013.

NIBSC/WHO (87/586) approximate value (U/mL) =  $0.008 \times \text{Human IFN-}\gamma \text{ value (pg/mL)}$ 

# **PLATE LAYOUT**



Note: Blank is a Blank Dilution Buffer hole.

10 / 12

**US and Canada**:

Tel: +1 800-810-0816

E-mail: order@acrobiosystems.com

Web: http://www.acrobiosystems.com

**Asia and Pacific**: **Tel**: +86 400-682-2521



# **TROUBLESHOOTING GUIDE**

Problem	Cause	Solution
Poor standard curve	* Inaccurate pipetting	* Check pipettes
Large CV * Inaccurate pipetting		* Check pipettes
	* Air bubbles in wells	* Remove bubbles in wells
High background	* Plate is insufficiently washed	* Review the manual for proper wash.
	* Contaminated wash buffer	* Make fresh wash buffer
Very low readings across	* Incorrect wavelengths	* Check filters/reader
the plate	* Insufficient development	* Increase development time
	time	
Samples are reading too * Samples contain cytokine * Di		* Dilute samples and run again
high, but standard curve levels above assay range		
looks fine		
Drift	* Interrupted assay set-up	* Assay set-up should be continuous - have all
	* Reagents not at room	standards and samples prepared appropriately
	temperature	before commencement of the assay
		* Ensure that all reagents are at room
		temperature before pipetting into the wells
		unless otherwise instructed in the antibody
		inserts