

T008-EN.02

# Anti-SARS-CoV-2 Antibody IgG Titer Serologic Assay Kit (Nucleocapsid Protein)

Pack Size: 96 tests

Catalog Number: RAS-T008

IMPORTANT: Please carefully read this manual before performing your experiment.

For Research Use Only. Not For Use In Diagnostic Or Therapeutic Procedu

HTTP://WWW.ACROBIOSYSTEMS.COM



#### **INTENDED USE**

This kit is developed for titer measurement of Anti-SARS-CoV-2 Antibody IgG (Nucleocapsid Protein) in human serum. It is intended for research use only (RUO).

#### PRINCIPLE OF THE ASSAY

The newly identified Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2) has posed a serious threat to human health. A rapid and effective Assay kit detecting the levels of Anti-SARS-CoV-2 in human serum can facilitate research on characterization of antibodies produced in response to SARS-CoV-2 infection.

This assay kit is used to measure the titer of Anti-SARS-CoV-2 Nucleocapsid Antibody IgG by employing an indirect ELISA. Immobilize SARS-CoV-2 Nucleocapsid Protein on the microplate. Then add the samples, incubate and wash the wells. Next add Secondary antibody HRP-Anti-Human IgG to the plate, incubate and wash the wells. Lastly load the substrate into the wells and monitor color development in proportion with the amount of antibody present. The reaction is stopped by the addition of a stop solution and the intensity of the absorbance can be measured at 450 nm. The OD Value reflects the amount of antibody bound.

#### MATERIALS PROVIDED

Catalog	Components	Size	Format	Storage	
		(96 tests)		Unopened	Opened
RAS008-C01	Pre-coated SARS-CoV-2 Nucleocapsid Protein Microplate	1 plate	Solid	2-8°C	2-8°C
RAS008-C02	SARS-CoV-2 Antibody Positive Control	100 µL	Liquid	2-8°C	2-8°C
RAS008-C03	SARS-CoV-2 Antibody Negative Control	100 µL	Liquid	2-8°C	2-8°C
RAS008-C04	HRP-Anti-Human IgG	200 µL	Liquid	2-8°C,avoid light	2-8°C,avoid light
RAS008-C05	10×Washing Buffer	50 mL	Liquid	2-8°C	2-8°C
RAS008-C06	Dilution Buffer	50 mL	Liquid	2-8°C	2-8°C
RAS008-C07	Substrate Solution	12 mL	Liquid	2-8°C,avoid light	2-8°C,avoid light

#### TABLE 1. MATERIALS PROVIDED

US and Canada:

Tel: +1 800-810-0816

Web: http://www.acrobiosystems.com

Asia and Pacific:

**Tel**: +86 400-682-2521

E-mail: order@acrobiosystems.com

2/5



RAS008-C08	Stop Solution	7 mL	Liquid	2-8°C	2-8°C	

## **REAGENTS/EQUIPMENT NEEDED BUT NOT SUPPLIED**

Single or dual wavelength microplate reader with 450 nm filter;

Centrifuge;

37 °C Incubator;

Single channel or multichannel pipettes with 10  $\mu$ L, 200  $\mu$ L and 1000  $\mu$ L precision;

10  $\mu$ L, 200  $\mu$ L and 1000  $\mu$ L pipette tips;

Test Tubes;

Graduated cylinder;

Deionized or distilled water for dilution;

### **STORAGE**

- 1. Unopened kit should be stored at 2°C-8°C upon receiving.
- 2. Find the expiration date on the outside packaging and do not use reagents past their expiration date.
- 3. The opened kit should be stored per components table. The shelf life is 30 days from the date of opening.

## **REAGENT PREPARATION**

Bring all reagents and samples to room temperature (20°C-25°C) before use. If crystals have formed in buffer solution,

warm until the crystals have completely dissolved and bring solution back to room temperature before use.

## **RECOMMENDED SAMPLE PREPARATION**

#### **1. Working fluid preparation**

1.1 Preparation of 1×Washing Buffer:

Dilute 50 mL 10×Washing Buffer with ultrapure water/deionized water to 500 mL.

1.2 Preparation of HRP-Anti-Human IgG working fluid:

Dilute HRP-Anti-Human IgG at 1:100 with Dilution Buffer. The prepared working fluid should avoid light. Please

prepare it for one-time use only.

1.3 Preparation of Positive Control and Negative Control working fluid and pre-treatment of samples:

Tel: +1 800-810-0816



T008-EN.02

It is recommended to dilute the samples, Positive Control and Negative Control from 1:400-1:12800 with Dilution Buffer.

#### 2. Add Samples

Add 100 µL diluted sample,Positive Control and Negative Control working fluid to the corresponding wells. Seal the plate with microplate sealing film and incubate at 37°C for 1.0 h.

#### 3. Washing

Remove the remaining solution by aspiration, add 300  $\mu$ L of 1×Washing Buffer to each well, gently tap the plate for 30 s, remove any remaining 1×Washing Buffer: by aspirating or decanting, invert the plate and blot it against paper towels. Repeat the wash step above for three times.

### 5.Add HRP-Anti-Human IgG

For all wells, add 100 µL HRP-Anti-Human IgG working solution.Seal the plate with microplate sealing film and incubate at 37°C for 1.0 h, avoid light.

#### 6. Washing

Repeat step 3.

#### 7. Substrate Reaction

Add 100 µL Substrate Solution to each well. Seal the plate with microplate sealing film and incubate at 37°C for 20 min, avoid light.

#### 8. Termination

Add 50 µL Stop Solution to each well, and tap the plate gently for 3 min to allow thorough mixing.

*Note:* the color in the wells should change from blue to yellow.

#### 9. Data Recording

Read the absorbance at 450 nm using UV/Vis microplate spectrophotometer.

*Note*: To reduce the background noise, substract the value read at  $OD_{450 \text{ nm}}$  with the value read at  $OD_{630 \text{ nm}}$ .

## **CUT-OFF VALUE IDENTIFICATION**

Cut-off value =0.1.

Note: The cut-off value can be determined by the end user.

#### **INTERPRETION OF RESULTS**

US and Canada: Asia and Pacific:

Tel: +1 800-810-0816 Tel: +86 400-682-2521



#### For determination of antibody titer:

Determination of antibody titer: the positive sample was diluted with a gradient, and the antibody titer of the sample corresponds to the highest dilution factor that still yields a positive reading.

## **LIMITATIONS OF THE PROCEDURE**

This kit is developed for detecting human serum of Anti-SARS-CoV-2 Antibody IgG (Nucleocapsid Protein). However, we do not have the LoQ (Limit of Quantitation) and ULMI (upper limit of measuring interval) and cutoff defined for semi-quantitative detection. Interested customer is recommended to establish the semi-quantitative detection procedure themselves.

## **PERFORMANCE**

Precision: Intra batch CV%<15%, Inter batch CV%<15%. Specificity: 93.8% (five samples show false positive ,n=80).

### **PRECAUTIONS**

- 1. This kit is for research use only and is not for use in diagnostic or therapeutic applications.
- 2. This kit should be used according to the provided instructions.
- 3. Do not mix reagents from different lots.
- 4. Bring all reagents and samples to room temperature (20°C-25°C) before use.
- If crystals have formed in the buffer solution, incubate until the crystals have completely dissolved. Before use, bring the solution back to room temperature.
- 6. This kit should be stored at 2°C-8°C.
- 7. Please prepare the working solution of each component according to the needs of the experiment, all prepared working solution is for one-time use and cannot be stored.

T008-EN.02

US and Canada: