

**Synonym**

TNFRSF17,CD269,BCM,BCMA

Source

Biotinylated Mouse BCMA, Fc,Avitag(BCA-M82F0) is expressed from human 293 cells (HEK293). It contains AA Met 1 - Thr 49 (Accession # [O88472-1](#)).

Predicted N-terminus: Met 1

Molecular Characterization

BCMA(Met 1 - Thr 49) O88472-1	Fc(Pro 100 - Lys 330) P01857	Avi
----------------------------------	---------------------------------	-----

This protein carries a human IgG1 Fc tag at the C-terminus, followed by an Avi tag (Avitag™).

The protein has a calculated MW of 34.4 kDa. The protein migrates as 38-43 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE) due to O-linked glycosylation which has been verified by O-Glycosidase and Neuraminidase digestion.

Labeling

Biotinylation of this product is performed using Avitag™ technology. Briefly, the single lysine residue in the Avitag is enzymatically labeled with biotin.

Protein Ratio

Passed as determined by the HABA assay / binding ELISA.

Endotoxin

Less than 1.0 EU per µg by the LAL method.

Purity

>90% as determined by SDS-PAGE.

Formulation

Lyophilized from 0.22 µm filtered solution in

Tris with Glycine, Arginine and NaCl, pH7.5 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.


Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.



This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

SDS-PAGE

Biotinylated Mouse BCMA, Fc,Avitag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90% (With [Star Ribbon Pre-stained Protein Marker](#)).


Bioactivity-ELISA**Discounts, Gifts,
and more!**




Immobilized Mouse BAFF, Mouse IgG2a Fc Tag (Cat. No. BAF-M5257) at 5 µg/mL (100 µL/well) can bind Biotinylated Mouse BCMA, Fc,Avitag (Cat. No. BCA-M82F0) with a linear range of 2-39 ng/mL (QC tested).

Immobilized Biotinylated Mouse BCMA, Fc,Avitag (Cat. No. BCA-M82F0) at 1 µg/mL (100 µL/well) on Recombinant Streptavidin (Cat. No. STN-N5116) precoated (0.5µg/well) plate, can bind Human APRIL, His,Flag Tag, active trimer (Cat. No. APL-H52D1) with a linear range of 10-78 ng/mL (Routinely tested).

Bioactivity-BLI



Loaded Biotinylated Mouse BCMA, Fc,Avitag (Cat. No. BCA-M82F0) on SA Biosensor, can bind Mouse BAFF, Mouse IgG2a Fc Tag, low endotoxin (Cat. No. BAF-M5257) with an affinity constant of 3.75 nM as determined in BLI assay (ForteBio Octet Red96e) (Routinely tested).



Loaded Biotinylated Mouse BCMA, Fc,Avitag (Cat. No. BCA-M82F0) on SA Biosensor, can bind Human APRIL, His,Flag Tag, active trimer (Cat. No. APL-H52D1) with an affinity constant of 2.19 nM as determined in BLI assay (ForteBio Octet Red96e) (Routinely tested).

Background

Tumor necrosis factor receptor superfamily member 17 (TNFRSF17) is also known as B-cell maturation protein (BCMA), CD antigen CD269, which is a member of the TNF-receptor superfamily. TNFRSF17 contains one TNFR-Cys repeat. TNFRSF17 is expressed in mature B-cells, but not in T-cells or monocytes. TNFRSF17 is receptor for TNFSF13B/BLyS/BAFF and TNFSF13/APRIL. TNFRSF17 promotes B-cell survival and plays a role in the regulation of humoral immunity. TNFRSF17 can activate NF-kappa-B and JNK.

Clinical and Translational Updates

