

Mouse Cadherin-17 / CDH17 Protein, Fc Tag (MALS verified)

Catalog # CA7-M5256



BIOSYSTEMS
Acro
Surprise Inside!

Synonym

Cadherin-17, CDH17, HPT-1, LI-cadherin

Source

Mouse Cadherin-17 Protein, Fc Tag(CA7-M5256) is expressed from human 293 cells (HEK293). It contains AA Phe 26 - Met 786 (Accession # [Q9R100-1](#)).

Predicted N-terminus: Phe 26

Molecular Characterization

Cadherin-17(Phe 26 - Met 786) Q9R100-1	Fc(Pro 100 - Lys 330) P01857
---	---------------------------------

This protein carries a human IgG1 Fc tag at the C-terminus.

The protein has a calculated MW of 110.8 kDa. The protein migrates as 50 kDa and 110-130 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE) due to glycosylation.

Endotoxin

Less than 1.0 EU per µg by the LAL method.

Purity

>90% as determined by SDS-PAGE.

>90% as determined by SEC-MALS.

Formulation

Lyophilized from 0.22 µm filtered solution in 20 mM Tris, 150 mM NaCl, pH8.0 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage


For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.


SDS-PAGE



Mouse Cadherin-17 Protein, Fc Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90% (With [Star Ribbon Pre-stained Protein Marker](#)).

Bioactivity-ELISA

SEC-MALS



The purity of Mouse Cadherin-17 Protein, Fc Tag (Cat. No. CA7-M5256) is more than 90% and the molecular weight of this protein is around 240-290 kDa verified by SEC-MALS.

[Report](#)

Discounts, Gifts,
and more!




➤ www.acrobiosystems.com

12/8/2023

Mouse Cadherin-17 / CDH17 Protein, Fc Tag (MALS verified)

Catalog # CA7-M5256



Immobilized Mouse Cadherin-17 Protein, Fc Tag (Cat. No. CA7-M5256) at 1 µg/mL (100 µL/well) can bind various dilution ratio of Anti-Cadherin 17 (QC tested).

Background

Cadherin-17, also known as liver-intestine (LI) Cadherin, belongs to the cadherin family of calcium-dependent cell adhesion molecules. In vivo studies showed CDH17 knockout resulted in apoptotic PC tumor death through activating caspase-3 activity. Taken together, CDH17 functions as an oncogenic molecule critical to PC growth by regulating tumor apoptosis signaling pathways and CDH17 could be targeted to develop an anti-PC therapeutic approach.

Clinical and Translational Updates

Please contact us via TechSupport@acrobiosystems.com if you have any question on this product.

Discounts, Gifts,
and more!

