



## Synonym

CD5,LEU1

## Source

Biotinylated Human CD5, His,Avitag, premium grade(CD5-H82E5) is expressed from human 293 cells (HEK293). It contains AA Arg 25 - Pro 372 (Accession # [P06127-1](#)).

Predicted N-terminus: Arg 25

*It is produced under our rigorous quality control system that incorporates a comprehensive set of tests including sterility and endotoxin tests. Product performance is carefully validated and tested for compatibility for cell culture use or any other applications in the early preclinical stage. When ready to transition into later clinical phases, we also offer a custom GMP protein service that tailors to your needs. We will work with you to customize and develop a GMP-grade product in accordance with your requests that also meets the requirements for raw and ancillary materials use in cell manufacturing of cell-based therapies.*

## Molecular Characterization

|                                   |          |     |
|-----------------------------------|----------|-----|
| CD5(Arg 25 - Pro 372)<br>P06127-1 | Poly-his | Avi |
|-----------------------------------|----------|-----|

This protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (Avitag™).

The protein has a calculated MW of 42.2 kDa. The protein migrates as 50-55 kDa under reducing (R) condition, and 45-50 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under non-reducing (NR) condition (SDS-PAGE) due to glycosylation.

## Labeling

*Biotinylation of this product is performed using Avitag™ technology. Briefly, the single lysine residue in the Avitag is enzymatically labeled with biotin.*

## Biotinylation

As determined by Quantative ELISA binding assay against streptavidin.

## Endotoxin

Less than 0.01 EU per µg by the LAL method.

## Host Cell Protein

<0.5 ng/µg of protein tested by ELISA.

## Host Cell DNA

<0.02 ng/µg of protein tested by qPCR.

## Sterility

Negative

## Mycoplasma

Negative.

## Purity

>95% as determined by SDS-PAGE.

## Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

## Reconstitution

Please see Certificate of Analysis for specific instructions.

*For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.*

## Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower.

*Please avoid repeated freeze-thaw cycles.*


This product is stable after storage at:

- -20°C to -70°C for 24 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

## SDS-PAGE



Discounts, Gifts,  
and more!






Biotinylated Human CD5, His,Avitag, premium grade on SDS-PAGE under reducing (R) and non-reducing (NR) conditions. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95% (With [Star Ribbon Pre-stained Protein Marker](#)).

**Bioactivity-ELISA**



Immobilized Biotinylated Human CD5, His,Avitag, premium grade (Cat. No. CD5-H82E5) at 1 µg/mL (100 µL/well) on streptavidin (Cat. No. STN-N5116) precoated (0.5 µg/well) plate can bind Mouse Anti-CD5 Antibody, Mouse IgG1 with a linear range of 0.1-1.6 ng/mL (QC tested).

Immobilized Biotinylated Human CD5, His,Avitag, premium grade (Cat. No. CD5-H82E5) at 1 µg/mL (100 µL/well) on streptavidin (Cat. No. STN-N5116) precoated (0.5 µg/well) plate can bind Monoclonal Anti-Human CD5 Antibody, Mouse IgG1 (2C8A7) with a linear range of 0.06-4 ng/mL (Routinely tested).



Immobilized Biotinylated Human CD5, His,Avitag, premium grade (Cat. No. CD5-H82E5) at 1 µg/mL (100 µL/well) on streptavidin (Cat. No. STN-N5116) precoated (0.5 µg/well) plate can bind Monoclonal Anti-Human CD5

Discounts, Gifts, and more!



# Biotinylated Human CD5 Protein, His,Avitag™, premium grade

Catalog # CD5-H82E5



BIOSYSTEMS  
**Acro**

Antibody, Mouse IgG1 (9F8F11) (Cat. No. CD5-MY2044) with a linear range of 0.06-4 ng/mL (Routinely tested).

## Background

T-cell surface glycoprotein CD5 is also known as Lymphocyte antigen T1/Leu-1 and LEU1, which is phosphorylated on tyrosine residues by LYN, so CD5 can create binding sites for PTPN6/SHP-1. CD5 may act as a receptor in regulating T-cell proliferation. CD5 is expressed at various developmental and activation stages on human B cells. CD5 is a well established negative regulator of TCR and BCR signalling. CD5-positive cells may also prevent the emergence of autoimmunity by provision of cytokines such as IL-10. Development, selection and function of different B- and T-cell subsets or their preferential survival may be directly or indirectly dependent on different glycan structures associated with CD5 or CD5-like molecules.

## Clinical and Translational Updates

Discounts, Gifts,  
and more!



➤ [www.acrobiosystems.com](http://www.acrobiosystems.com)

8/12/2024