



## Synonym

AITR,GITR,TNFRSF18,CD357

## Source

Human GITR, Fc Tag(GIR-H5254) is expressed from human 293 cells (HEK293). It contains AA Gln 26 - Glu 161 (Accession # [Q9Y5U5-1](#)).

Predicted N-terminus: Gln 26

## Molecular Characterization

GITR(Gln 26 - Glu 161) Q9Y5U5-1	Fc(Pro 100 - Lys 330) P01857
------------------------------------	---------------------------------

This protein carries a human IgG1 Fc tag at the C-terminus.

The protein has a calculated MW of 41.0 kDa. The protein migrates as 45-47 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

## Endotoxin

Less than 1.0 EU per µg by the LAL method.

## Purity

>95% as determined by SDS-PAGE.

## Formulation

Lyophilized from 0.22 µm filtered solution in

Tris with Glycine, Arginine and NaCl, pH7.5 with trehalose as protectant.

Contact us for customized product form or formulation.

## Reconstitution

Please see Certificate of Analysis for specific instructions.

*For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.*

## Storage


For long term storage, the product should be stored at lyophilized state at -20°C or lower.

*Please avoid repeated freeze-thaw cycles.*

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

## SDS-PAGE



Human GITR, Fc Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%.


## Bioactivity-ELISA

Discounts, Gifts,  
and more!





Surprise Inside!



Immobilized Human GITR, Fc Tag (Cat. No. GIR-H5254) at 2 µg/mL (100 µL/well) can bind Biotinylated Human GITR Ligand, Avitag, Fc Tag (Cat. No. GIL-H82F8) with a linear range of 4-63 ng/mL (QC tested).

## Background

Glucocorticoid-induced TNFR-related protein (GITR) is also known as Tumor necrosis factor receptor superfamily member 18 (TNFRSF18), activation-inducible TNFR family receptor (AITR), CD antigen CD357, which is a member of the tumor necrosis factor receptor (TNF-R) superfamily. GITR is receptor for TNFSF18, which seems to be involved in interactions between activated T-lymphocytes and endothelial cells and in the regulation of T-cell receptor-mediated cell death. GITR also mediated NF-kappa-B activation via the TRAF2/NIK pathway.

## Clinical and Translational Updates

Discounts, Gifts,  
and more!

