



Source

Monoclonal Anti-IFN γ Antibody, Mouse IgG2a (8C5F8) is a Mouse monoclonal antibody produced from a hybridoma created by fusing SP2/0 myeloma and Mouse B-lymphocytes.

Clone

8C5F8

Species

Mouse

Isotype

Mouse IgG2a/kappa

Conjugate

Unconjugated

Antibody Type

Hybridoma Monoclonal

Reactivity

Human

Immunogen

Recombinant Human IFN γ is expressed from human HEK293 cells.

Specificity

This product is a specific antibody specifically reacts with IFN- γ , Human. No cross-reactivity is detected with other human cytokines, including IL-2, IL-4, IL-6, IL-10, GM-CSF and TNF-alpha.

Application

Application	Recommended Usage
ELISA	0.2-25 ng/mL

Purity

>95% as determined by SDS-PAGE.

Purification

Protein A purified/ Protein G purified

Formulation

Lyophilized from 0.22 μ m filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage

For long term storage, the product should be stored at lyophilized state at -20 $^{\circ}$ C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20 $^{\circ}$ C to -70 $^{\circ}$ C for 12 months in lyophilized state;
- -70 $^{\circ}$ C for 3 months under sterile conditions after reconstitution.


SDS-PAGE

Discounts, Gifts,
and more!



Monoclonal Anti-IFN γ Antibody, Mouse IgG2a (8C5F8)

Catalog # IFN-S120




Monoclonal Anti-IFN γ Antibody, Mouse IgG2a (8C5F8) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%.

Bioactivity-ELISA


Monoclonal Anti-IFN γ Antibody, Mouse IgG2a (8C5F8) ELISA

0.2 μ g of Human IFN-gamma / IFNG Protein, premium grade per well



Immobilized Human IFN-gamma, premium grade (Cat. No. IFG-H4211) can bind Monoclonal Anti-IFN γ antibody, Mouse IgG2a (8C5F8) (Cat. No. IFN-S120) with a linear range of 0.02-25 ng/mL (QC tested).

Bioactivity-SPR



Monoclonal Anti-IFN γ antibody, Mouse IgG2a (8C5F8) (Cat. No. IFN-S120) captured on CM5 chip via anti-mouse antibodies surface can bind Human IFN-gamma, premium grade (Cat. No. IFG-H4211) with an affinity constant of 0.933 nM as determined in a SPR assay (Biacore 8K) (Routinely tested).

Discounts, Gifts,
and more!



Monoclonal Anti-IFN γ Antibody, Mouse IgG2a (8C5F8)

Catalog # IFN-S120



Background

Interferon-gamma (IFN- γ /IFNG) is a dimerized soluble cytokine that is the only member of the type II class of interferon. This interferon was originally called macrophage-activating factor, a term now used to describe a larger family of proteins to which IFN- γ belongs. IFN-gamma has been used in a wide variety of clinical indications. Interferon-gamma (IFN γ) is a central regulator of the immune response and signals via the Janus Activated Kinase (JAK)-Signal Transducer and Activator of Transcription (STAT) pathway. Interferon gamma has broader roles in activation of innate and adaptive immune responses to viruses and tumors, in part through upregulating transcription of genes involved in cell cycle regulation, apoptosis, and antigen processing/presentation. Despite this, rodent and human trophoblast cells show dampened responses to IFNG that reflect the resistance of these cells to IFNG-mediated activation of major histocompatibility complex (MHC) class II transplantation antigen expression.

Clinical and Translational Updates

Discounts, Gifts,
and more!

