

Anti-SARS-CoV-2 Spike RBD Neutralizing Antibody, Chimeric mAb, Human IgG1 (AM122) (Trehalose free)

Catalog # S1N-M12A1



Source

Anti-SARS-CoV-2 Spike RBD Antibody, Chimeric mAb, Human IgG1 (AM122) is a chimeric monoclonal antibody recombinantly expressed from HEK293 cells, which combines the variable region of a mouse monoclonal antibody with human IgG1 constant domain. The mouse monoclonal antibody was obtained from a mouse immunized with recombinant SARS-CoV-2 Spike S1 protein. *As verified in competitive ELISA-based and pseudovirus-based neutralization assay, this chimeric monoclonal can potently neutralize all SARS-CoV-2 Variants of Concern (VOCs), including Alpha (B.1.1.7), Beta (B.1.351), Gamma (P.1) and Delta (B.1.617.2).*

Clone

AM122

Isotype

Human IgG1 | Human Kappa

Conjugate

Unconjugated

Antibody Type

Recombinant Monoclonal

Reactivity

Virus

Specificity

This product is a specific antibody against SARS-CoV-2 Spike protein RBD domain. No cross-reactivity is detected with Spike protein RBD domain of other coronaviruses, including SARS-CoV, MERS-CoV, HCoV-229E, HCoV-NL63, HCoV-OC43 and HCoV-HKU1.

Application

| Application | Recommended Usage |
|-------------|-------------------|
| ELISA | 0.1-2.5 ng/mL |

Purity

>95% as determined by SDS-PAGE.

>90% as determined by SEC-MALS.

Purification

Protein A purified/ Protein G purified

Formulation

Supplied as 0.2 µm filtered solution in PBS, pH7.4 .

Contact us for customized product form or formulation.

Shipping

This product is supplied and shipped with dry ice, please inquire the shipping cost.

Storage

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- For long term storage, the product is stable for up to 3 years at -70°C from date of receipt;
- For short term storage, the product is stable for up to 12 months at 2-8°C from date of receipt.

SDS-PAGE


SEC-MALS

Discounts, Gifts,
and more!




Anti-SARS-CoV-2 Spike RBD Neutralizing Antibody, Chimeric mAb, Human IgG1 (AM122) (Trehalose free)

Catalog # S1N-M12A1




Anti-SARS-CoV-2 Spike RBD Antibody, Chimeric mAb, Human IgG1 (AM122) on SDS-PAGE under reducing (R) and non-reducing (NR) conditions. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%.


Bioactivity-ELISA



Immobilized SARS-CoV-2 S1 protein, His Tag (Cat. No. S1N-C52H2) at 2 µg/mL (100 µL/well) can bind Anti-SARS-CoV-2 Spike RBD Antibody, Chimeric mAb, Human IgG1 (AM122) (Cat. No. S1N-M12A1) with a linear range of 0.08-2.5 ng/mL (QC tested).




Anti-SARS-CoV-2 Spike RBD Neutralizing antibody (Cat.No. S1N-M12A1) neutralizes SARS-CoV-2 Spike RBD by inhibiting RBD: ACE2 interaction. The ACE2-coated plate is incubated with the wild type (WT) RBD or B.1.1.7, B.1.351, P.1, B.1.617.1, B.1.617.2 mutant and treated with the neutralizing antibody at increasing concentration. Percent inhibition is calculated based on the OD value.




The purity of Anti-SARS-CoV-2 Spike RBD Antibody, Chimeric mAb, Human IgG1 (AM122) (Cat. No. S1N-M12A1) is more than 90% and the molecular weight of this protein is around 135-150 kDa verified by SEC-MALS.

[Report](#)



Detection SARS-CoV-2 S protein RBD by Sandwich ELISA Assay.

Immobilized Anti-SARS-CoV-2 Spike S1 Antibody, chimeric mAb (Cat. No. S1N-M12A1) at 2 µg/mL (100 µL/well) can bind S protein RBD. And then add Anti-SARS-CoV-2 Spike S1 Antibody (Cat. No. S1N-M13A1) (Biotinylated in house) at 1:5000. Detection was performed using HRP-conjugated streptavidin with sensitivity of 25 pg/mL (Routinely tested).



Detection of HRP-SARS-CoV-2 Spike RBD(B.1.1.529/Omicron), His Tag titer by competitive-ELISA Assay. Serial dilutions of Anti-SARS-CoV-2 Spike RBD Neutralizing Antibody, Chimeric mAb, Human IgG1 (AM122) (Trehalose free) (Cat. No. S1N-M12A1) were added into Human ACE2 / ACEH Protein, Fc Tag (Cat. No. AC2-H5257): HRP-SARS-CoV-2 Spike RBD(B.1.1.529/Omicron), His Tag binding reactions. The half maximal inhibitory concentration (IC50) is 0.340 µg/mL (Routinely tested).

Discounts, Gifts,
and more!



Anti-SARS-CoV-2 Spike RBD Neutralizing Antibody, Chimeric mAb, Human IgG1 (AM122) (Trehalose free)

Catalog # S1N-M12A1



Background

It's been reported that Coronavirus can infect the human respiratory epithelial cells through interaction with the human ACE2 receptor. The spike protein is a large type I transmembrane protein containing two subunits, S1 and S2. S1 mainly contains a receptor binding domain (RBD), which is responsible for recognizing the cell surface receptor. S2 contains basic elements needed for the membrane fusion. The S protein plays key parts in the induction of neutralizing-antibody and T-cell responses, as well as protective immunity.

Clinical and Translational Updates

Discounts, Gifts,
and more!

