

Synonym

Spike,S protein RBD,Spike glycoprotein Receptor-binding domain,S glycoprotein RBD,Spike protein RBD

Source

SARS-CoV-2 Spike RBD Protein, His Tag (BA.2.86/V483del/USA) (SPD-C5247) is expressed from human 293 cells (HEK293). It contains AA Arg 319 -Lys 537 (Accession # QHD43416.1

(I332V,G339H,K356T,S371F,S373P,S375F,T376A,R403K,D405N,R408S,K417 N,N440K,V445H,G446S,N450D,L452W,N460K,S477N,T478K,N481K,V483de l,E484K,F486P,Q498R,N501Y,Y505H)). The spike mutations are identified on the SARS-CoV-2 Omicron variant (Pango lineage: BA.2.86/V483del/USA). Predicted N-terminus: Arg 319

Molecular Characterization

This protein carries a polyhistidine tag at the C-terminus.

The protein has a calculated MW of 26.7 kDa. The protein migrates as 38-43 kDa when calibrated against <u>Star Ribbon Pre-stained Protein Marker</u> under reducing (R) condition (SDS-PAGE) due to glycosylation.

Endotoxin

Less than $1.0\ EU$ per μg by the LAL method.

Purity

>90% as determined by SDS-PAGE.

Formulation

Lyophilized from $0.22~\mu m$ filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage


For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.


SDS-PAGE



SARS-CoV-2 Spike RBD Protein, His Tag (BA.2.86/V483del/USA) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90% (With <u>Star Ribbon Pre-stained Protein Marker</u>).

Bioactivity-ELISA

SEC-MALS




The purity of SARS-CoV-2 Spike RBD Protein, His Tag (BA.2.86/V483del/USA) (Cat. No. SPD-C5247) is more than 85% and the molecular weight of this protein is around 35-45 kDa verified by SEC-MALS. Report

SARS-CoV-2 Spike RBD Protein, His Tag (BA.2.86/V483del/USA) (MALS verified)

Catalog # SPD-C5247



SARS-CoV-2 Spike RBD Protein, His Tag (BA.2.86/V483del/USA) ELISA 0.5 µg of Human ACE2, Fc Tag per well




Immobilized Human ACE2, Fc Tag (Cat. No. AC2-H5257) at 5 μ g/mL (100 μ L/well) can bind SARS-CoV-2 Spike RBD Protein, His Tag

(BA.2.86/V483del/USA) (Cat. No. SPD-C5247) with a linear range of 1.953-

31.25 ng/mL (QC tested).

SARS-CoV-2 Spike RBD Protein, His Tag (BA.2.86/V483del/USA) ELISA 0.1 µg of SARS-CoV-2 Spike RBD Protein, His Tag (BA.2.86/V483del/USA) per well



Immobilized SARS-CoV-2 Spike RBD Protein, His Tag (BA.2.86/V483del/USA) (Cat. No. SPD-C5247) at 1 μ g/mL (100 μ L/well) can bind Human ACE2, Fc Tag (Cat. No. AC2-H5257) with a linear range of 0.122-3.906 ng/mL (Routinely tested).

Background

It's been reported that coronavirus can infect the human respiratory epithelial cells through interaction with the human ACE2 receptor. The spike protein is a large type I transmembrane protein containing two subunits, S1 and S2. S1 mainly contains a receptor binding domain (RBD), which is responsible for recognizing the cell surface receptor. S2 contains basic elements needed for the membrane fusion. The S protein plays key parts in the induction of neutralizing-antibody and T-cell responses, as well as protective immunity.

Clinical and Translational Updates

Please contact us via <u>TechSupport@acrobiosystems.com</u> if you have any question on this product.